

Contrast RED Solution

for Immunohistochemistry



Catalog No.
71-00-05

Size
50 mL

DESCRIPTION

Contrast RED Solution for Immunohistochemistry contains a nuclear fast red stain. It is designed for use as a counterstain with the HistoMark immunohistochemical staining systems.

FORM/STORAGE/STABILITY

Contrast RED is a 50 mL solution. Store tightly capped, refrigerated at 2 - 8°C. Stable for at least one year stored at 2 - 8°C.

USE

Product is provided at working concentration; no dilution is required.

After immunoassay is complete and tissue sections or other cellular preparations have been reacted with substrate solution, rinse slide with running tap water for 2 - 3 minutes. Counterstain specimen with Contrast RED Solution for 5 - 10 minutes.

After incubation with Contrast RED, rinse specimens thoroughly in reagent quality water. Slides may be dehydrated through alcohol to xylene or air dried and mounted in aqueous or xylene-based mounting media. **DO NOT USE XYLENE OR ALCOHOL WITH THE HISTOMARK BLUE STAINING SYSTEM!**

IMMUNOASSAY KIT REAGENTS

Contrast RED Solution for Immunohistochemistry is a component of the HistoMark® BLUE kit for Alkaline Phosphatase Conjugate Reporter Systems.

HistoMark RED	Catalog Number 55-69-00
HistoMark BLUE	Catalog Number 55-70-00

RELATED PRODUCTS

HistoMark® RED Phosphatase System	Cat.No.55-69-00
HistoMark® BLUE Phosphatase System	Cat.No.55-70-00
HistoMark® BLACK Peroxidase System	Cat.No.54-75-00
HistoMark® ORANGE Peroxidase System	Cat.No.54-74-00

See KPL catalog for a complete list of antibodies, substrates and complete systems for microwell ELISA, membrane blotting and immunohistochemical applications.

L-535-02

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product. No claims beyond replacement of unacceptable material or refund of purchase price shall be allowed.

Gaithersburg, MD • USA

301-948-7755 • 800-638-3167 • FAX 301-948-0169 • www.kpl.com